

**FDA INSERT**  
**Important Information Regarding**  
**BU-1MI**

**Intended Use:**

Blood/urine specimen collection and transport

**Contents:**

FSD-93 blood/urine alcohol report form  
blood tube protector insert  
10 ml blood collection tubes containing 100 mg of sodium fluoride  
and 20 mg of potassium oxalate  
alcohol-free prep pad  
urine specimen collection bottles  
ziplock bag containing liquid absorbing sheet and important note  
specimen labels  
kit shipping seal

**Warnings and Precautions:**

Urine and blood samples should be handled and processed  
as if they are potentially infectious.

**Waste Disposal Instructions:**

- **Needles:**  
Dispose of any needles in an approved Sharps container.
- **Blood Tubes:**  
Dispose of all tubes using safe laboratory procedures as  
outlined in bio-safety and microbiological and bio-medical  
laboratories as HH Publication CDC 84-8395.

This product information sheet is included to comply  
with FDA Regulations.

BU1MI: FDA.2 8/97

**IMPORTANT!!**

**Place blood tube holder containing  
filled blood tubes along with filled  
specimen bottles inside ziplock bag,  
then squeeze out all excess air and  
close bag.**

***NOTE: DO NOT REMOVE LIQUID  
ABSORBING SHEET.***

BU1MI: IN.1 3/96